

CEREC 4.0x cheat sheet

General

Ctrl + N	New case
Ctrl + S	Save
Ctrl + Alt + S	Save as
Ctrl + I	Import
Ctrl + E	Export
F1	Help
Ctrl + F1	Administration phase
Ctrl + F2	Acquisition phase
Ctrl + F3	Model phase
Ctrl + F4	Design phase
Ctrl + F5	Mill phase
Ctrl + F6	Configuration options
Alt + F4	Exit program

Acquisition phase

1	Lower jaw
2	Upper jaw
3	Buccal
4	BioRef lower
5	BioRef upper
6	BioCopy lower
7	BioCopy upper
Up arrow	Go to the next image catalog to the right
Down arrow	Go to the next image catalog to the left
Space	Toggle camera on/off

Design phase

Space	(in Form tool) Toggle between Remove, Smooth, Add
Space	(in Move tool) Toggle between Position/Rotate
Space	Toggle between Anatomical and Circular shape tool

Mill phase

Space	(in Move tool) Toggle between Position/Rotate tool
-------	--

Source: pattersonsupport.custhelp.com

Last modification: 6/17/2022 8:20:14 AM

More information: defkey.com/cerec-4-0x-shortcuts

[Customize this PDF...](#)