

Merge PACS Workstation cheat sheet

Non-image review actions

Alt + A	System-wide: Launch Merge Messenger
Alt + K	System-wide: Display assigned keyboard shortcut
Alt + L	Open new Local Study
Alt + M	Open Management Pages
Alt + Q	Open new Query Search
Alt + R	Open new RealTime Worklist
Alt + T	Open new Teaching Worklist
Alt + V	Open new Recently Viewed
Ctrl + L	Switch to Local Study
Ctrl + P	Switch to Patient Record
Ctrl + Q	Switch to Query Search
Ctrl + R	Switch to RealTime Worklist
Ctrl + T	Switch to Teaching Worklist
Ctrl + V	Switch to Recently Viewed
Ctrl + X	System-Wide: Close application

Last modification: 8/16/2020 8:28:27 AM

More information: defkey.com/merge-pacs-workstation-shortcuts

[Customize this PDF...](#)