



Affinity 2.0 & Compact cheat sheet

Tone

Up arrow	Channel 1 down
Down arrow	Channel 1 up
Left arrow	Frequency down
Right arrow	Frequency up
Page Up	Channel 2 down
Page Down	Channel 2 up
Shift + A	Single/dual audiogram
Shift + H	HF on/off
Shift + Z	HFz on/off
Shift + H	MF on/off

Speech graph

Up arrow	Correct
Down arrow	Incorrect
Left arrow	Channel 1 down
Right arrow	Channel 1 up
Page Up	Channel 2 up
Page Down	Channel 2 down

Speech table

Up arrow	Channel 1 up
Down arrow	Channel 1 down
Left arrow	Incorrect
Right arrow	Correct
Page Up	Channel 2 up
Page Down	Channel 2 down

Save & Print

Alt + S	Save session
Alt + X	Save & Exit
Alt + P	Print

General

F1	TF
F2	Tone test
F3	Speech test
F4	Weber test

F5	MHA test
F6	HLS test
F7	MLD test
F8	Knaster test
F9	Q'sin test
F10	SISI test
F11	Monitor on/off
F12	TB
5	Ch1 man/Rev
6	♪/♪♪
8	Ch2 man/Rev
9	Speech
9	Play/Pause
0	Speech
0	Stop
W	Ch2 on/off
E	+20 dB
R	Ch1 right
L	Ch1 left
D	0
F	1
G	2
H	3
J	4
C	CNT
B	Bone
N	No resp.
Space or Left Ctrl	Channel 1 stimulate
Right Ctrl	Channel 2 stimulate
Enter	Store heard
Del	Delete threshold

Source: interacoustics.com

Last modification: 02.03.2023 11:22:31

More information: defkey.com/pl/affinity-2-0-compact-hearing-aid-shortcuts

[Customize this PDF...](#)